




CMA[®]

CORONAL • MEDIAN • APICAL

SAFETY
and **SIMPLICITY**



CMA: a revolutionary system for your everyday practice

■ A reduced number of instruments.

4 nickel-titanium rotary endodontic instruments for a simplified sequence.

■ Only one sequence for root canal shaping and removal of filling material.

Refill of 6 NiTi rotary instruments			Taper	Tip ø	Ring	Total length
C	Coronal		10%	0.25 mm	White	17 mm
M	Median		6%	0.25 mm	Yellow	21 mm 25 mm
A	Apical 1		4%	0.20 mm	Red	21 mm 25 mm 29 mm
	Apical 2		6%	0.20 mm	Blue	21 mm 25 mm 29 mm

■ Additional instruments.

Refill of 6 NiTi rotary instruments		Taper	Tip ø	Ring	Total length
Presequence		2%	0.17 mm	Mauve	21 mm 25 mm
Apical 3		6%	0.30 mm	Green	21 mm 25 mm 29 mm
Apical 4		7%	0.40 mm	Black	21 mm 25 mm 29 mm

■ Only one sequence for simple, curved and narrow canals.

Clinical case by Dr. Stéphane Simon



(1) Pre-op X-ray



(2)



(2) (3) Post-op X-ray

Clinical case by Prof. Roger Rebeiz




(1) Maxillary first molar: narrow canals; dentine heavily mineralised.



(2) Constant and regular conicity obtained using CMA instruments.

Clinical case by Dr. Roland Arsan



(1) Pre-op X-ray: maxillary premolars and molars.



(2) The same CMA sequence was used for all 11 canals. All root canal trajectories respected.



(3) Same teeth, seen from another angle.
Only 4 instruments were needed to shape all these canals.


■ **Safer,
more secure
and reliable
rotary
instruments.**




>> Cross-section of 3 cutting edges ensures **greater cutting efficiency.**




>> Non-cutting tip: **ensures the trajectory of the root canal is respected.**



>> Their shape facilitates the **removal of debris from the canal.**




>> Short handle: **improves the access to molars.**



■ **Guaranteed Safety**

>> Stronger hard wearing NiTi alloy.
>> Deformation immediately visible to the naked eye.




Operating instructions for root canal shaping

Pre-op X-ray.



Opening
of the
pulp chamber.




Step 1

Exploring the coronal 2/3 of the canal and making it permeable:

Instruments used:

Manual N° 10 and 15 steel files
and/or **Presequence**.

Used until they move freely in the canal.



Objective:

- Secure and prepare the coronal 2/3 of the canal for the use of CMA.



Step 2

Flaring the coronal and middle part of the canal:

Instruments used:

CORONAL

Flares out the coronal part of the canal.

MEDIAN

Enlarges the middle part of the canal.



Objective:

- Straighten and enlarge the root canal entries to provide continuity between the pulp chamber and the canals, so the tools can freely access the apical one-third.



Recommendations

- CMA instruments are recommended for use with **endo motors featuring an auto-reverse function**.
- Use NiTi in a portion of the canal which has been previously explored and prepared with a manual N° 15 file and/or Presequence.
- **Examine the instruments** before and after each use. Discard the tool if there is even the slightest deformation.
- **Time:** 5 to 10 seconds per rotary instruments.
- **Movement:** progression towards the apex by continuous short (1 to 3 mm) and rapid up-and-down strokes, finishing off with a sweeping movement up against the root canal walls.
- NiTi instruments **must never be forced**.
- **Clean** the instrument properly **after each removal** with a compress.
- The canal must be **copiously irrigated** with sodium hypochlorite each time the instrument has been inserted.
- Use of chelating gel is advised in order to ease work with the tools.
- If progress with **A1** is blocked, use the K 15 file and **M** again. If it is blocked when using **A2**, use the K 15 file and **A1** again.

Additional recommendations for retreatment

- NiTi instruments can remove fillings of materials which **can be softened with solvents**. They cannot be used to remove fillings of **insoluble resin paste**.
- As the tool moves further towards the apex, use less solvent and irrigate more liberally.


Step 3

Determining the working length:

Instruments used:

Manual steel files N° 10 and 15, and/or Presequence.

Used until they move freely in the canal.



Objectives:

- Determining the working length.
- Preparing for the nickel-titanium rotary instruments to pass safely all the way to the apical one-third.



After the use of the **C** and the **M**, interference in the cervical area and coronal curves is eliminated, which further frees the way to the apical one-third.

Step 4

Enlarging the apical portion of the canal:

Instruments used:

APICAL 1 is used on the whole length of the root canal, followed by **APICAL 2**.



Objectives:

- Preparing the apical portion.
- Making a space where the irrigation solutions can collect.
- Providing the apical taper which will give the best obturation.
- Keeping apical diameter small.



Prof. Roger Rebeiz

(1) Gutta points in place.

(2) (3) Root canal obturation seen from different angles.

Treatment

K-FILES



■ Exploring

CORONAL



■ Coronal flaring

MEDIAN



■ Enlarging the middle part

K-FILES



■ Working length

APICAL 1



■ Apical preparation

APICAL 2



■ Finishing before filling

Speed: 300 to 400 rpm – Torque: 1.3 to 2.3 N/cm

Operating instructions for retreatment



Preoperative X-ray

Access cavity.
This must provide a clear view of the root canal entries and adequate access.

- Cleaning out all traces of filling material.

- Use of an ultrasonic scaler is recommended.

- Application of an appropriate solvent in the pulp chamber.



Coronal Step

Removal of filling material and flaring of the coronal portion of the canal:

- **Manual penetration** (e.g. using a N° 10 steel file shortened by a few mm), to pierce the filling material and create a 2-3 mm channel.
- **Use of CORONAL** to widen the canal entries and remove filling from 2-3 mm with a withdrawal movement, pressing against the wall(s).
- **Irrigation** and solvent renewal.
- **Manual penetration** using a N° 15 steel file going a few millimetres deeper.
- **Use of MEDIAN** which works deeper than the CORONAL.
This tool works using traction.
- **Irrigation** and solvent renewal.

Apical Step

Preparation and filling removal the apical portion of the canal:


- A **pre-curved manual N° 15 file** is used to explore this portion of the canal.
If feasible, the length of the canal may be measured at this stage.
- **Copious irrigation.**
- **APICAL 1** is used on the whole length of the root canal that has been made accessible using the N° 15 file.
- **Copious irrigation.**
- **APICAL 2** is used to remove the filling material and clear the canal.
- **Check passage** to the apex using a N° 10 K file just beyond the extent of the work.

■ **Same instruments used for the removal of a root canal filling and for shaping.**

Clinical case by Prof. Roger Rebeiz



Clinical case by Dr. Roland Arsan



(1) Pre-op X-ray
(2) Post-op X-ray

CORONAL




■ Flaring and removal of filling material

MEDIAN



■ Flaring the middle part

APICAL 1



■ Removal of filling material and apical preparation

APICAL 2



■ Flaring and finishing

Speed: 400 to 600 rpm for C and M - 300 to 400 rpm for A1 and A2 – **Torque:** 3 N/cm for C and M - 1.3 to 2.3 N/cm for A1 and A2

Additional instruments and accessories

PRESEQUENCE



- **PRESEQUENCE** is used to explore the root canal or to determine the working length.
- The use of PRESEQUENCE improves the root canal vacuity.




Dr. Roland Arsan

- (1) PRESEQUENCE instrument in place
- (2) Root canal obturation


APICAL 3 and APICAL 4



- **APICAL 3** is used after APICAL 2 in medium and large root canals, or in the step-back technique.
- **APICAL 4** is used after APICAL 3 in large root canals, or in the step-back technique.



- (1) Pre-op X-ray
- (2) Post-op X-ray (Use of A3 and A4)



- (1) Pre-op X-ray
- (2) PRESEQUENCE allows you to reach the apex.
Double curvature treated with APICAL 3 with step-back technique.

Prof. Roger Rebelz




Paper Points



Gutta Percha Points



Endo-Clean Discs




Sequencer





Endo-Clean Holder

4 instruments, only one sequence!

Product	Content
CMA NiTi START KIT A Length 25 mm	• 4 NITI instruments: 1 CORONAL, 1 MEDIAN, 1 APICAL 1, 1 APICAL 2 • 2 manual steel files n° 10 & n° 15
CMA NiTi START KIT B Length 21 mm	• 4 NITI instruments: 1 CORONAL, 1 MEDIAN, 1 APICAL 1, 1 APICAL 2 • 2 manual steel files n° 10 & n° 15
CMA NiTi START KIT C Length 29 mm	• 4 NITI instruments: 1 CORONAL, 1 MEDIAN, 1 APICAL 1, 1 APICAL 2 • 2 manual steel files n° 10 & n° 15

Refill of 6 NiTi rotary instruments			Taper	Tip ø	Ring	Total length
C	Coronal		10%	0.25 mm	White	17 mm
M	Median		6%	0.25 mm	Yellow	21 mm 25 mm
A	Apical 1		4%	0.20 mm	Red	21 mm 25 mm 29 mm
	Apical 2		6%	0.20 mm	Blue	21 mm 25 mm 29 mm

Refill of 6 NiTi rotary instruments			Taper	Tip ø	Ring	Total length
Presequence			2%	0.17 mm	Mauve	21 mm 25 mm
Apical 3			6%	0.30 mm	Green	21 mm 25 mm 29 mm
Apical 4			7%	0.40 mm	Black	21 mm 25 mm 29 mm

Instruments developed by Prof. Roger Rebeiz and Dr. Roland Arsan.
Illustrations and comments by Prof. Youssef Haikel.